

## CD19 (4G7)

Type	Size	Catalog number
Unconjugated	100µg	102901
	500µg	102903
FITC	25 tests	102914
	100 tests	102915
	200 tests	102916
PE	25 tests	102924
	100 tests	102925
	200 tests	102926
APC	25 tests	102944
	100 tests	102945
	200 tests	102946
PerCP	25 tests	102934
	100 tests	102935
	200 tests	102936
PerCP-Cyanine 5.5	25 tests	102964
	100 tests	102965
	200 tests	102966
iFluor™ 700	25 tests	1029194
	100 tests	1029195
	200 tests	1029196
Biotin	100µg	102951

**Antigen:** CD19  
**Immunogen:** Human CCL (chronic lymphocytic leukemia) cells  
**Host/Isotype:** Mouse, IgG1,κ  
**Reactivity:** Human  
**Purity:** >90% pure tested via polyacrylamide gel electrophoresis (PAGE)  
**Formulation:** PBS, pH7.2, 0.09%NaN<sub>3</sub> (unconjugated, Biotin)  
 PBS, pH7.2, 0.09% NaN<sub>3</sub> and 0.2% (w/v) BSA (conjugated)  
**Storage:** Store at 2-8°C and protected from prolonged exposure to light. **Do not freeze.**  
**Applications:** Flow Cytometry

❖ iFluor and mFluor are trademarks of AAT Bioquest, Inc.

## Application Information

Each lot of these antibodies has been pre-titrated and tested by flow cytometric analysis of human PBMCs such that 0.5µg (unconjugated, Biotin) or 5µl (conjugated) of these products are sufficient for staining 1 million cells in a 100µl staining volume or 100µl of whole blood. It is recommended to titrate antibody reactivity empirically for optimal performance.

**Antigen Information**

The clone 4G7 recognizes a 90-kDa CD19 antigen that is present on human B lymphocytes. The CD19 antigen is present on approximately 7 to 23% of human peripheral blood lymphocytes at all stages of B cell maturation but is lost on terminally differentiated plasma cells. CD19 does not react with resting or activated T lymphocytes, granulocytes, or monocytes.

**References**

1. Kugler M, et al. 2009. Protein Eng Des Sel. 22:135.
2. Lemmers B, et al. 2000. Leukemia. 14:2103.
3. Meeker TC, et al. 1984. Hybridoma. 3:305.
4. Strickler JG, et al. 1988. Cancer. 61:1782.

**Terms and Conditions**

This product is for research use only (RUO) and not intended for diagnostic testing.